

Updates and Reminders

Edward Livant, Wayne Hall, and Lorna Rabe
MTN Core Laboratory
Pittsburgh, PA





VOICE shipping and close out

- Thanks to the sites for all the efforts to get out shipments quickly to expedite endpoint verification for the DSMB.
- Please let the NL know if any documents are needed to be sent with shipping requests to expedite the process.
- Please carefully review instructions in the shipping emails



VOICE shipping and close out

- As study closeout begins, there will be pressure on the sites and NL to verify endpoints and perform HIV testing QA
- Once study visits are complete, other specimens such as cervical and vaginal swabs will be shipped.
- Gram stain shipments will be ongoing...



Sample Destruction Procedure

- Specimens collected from screening participants who did not enroll.
 - New LDMS report
 - Cross check with your site screening and enrollment report
 - Must be done before shipping specimens to the NL



Sample Destruction Procedure

- Specimens collected from participants who enrolled but did not consent to long term storage
 - Do not destroy these samples until notified by SCHARP and the NL
 - SCHARP will generate a list of PTIDS to be destroyed



EQA updates

- hCG
 - Switching from CAP to Accutest in 2012
 - Accutest sends 15 samples per year.
 - Some samples have low quantities of hCG
- Syphilis EQA vs. study samples
 - Report both the RPR and TPHA
- Save CAP and Accutest samples
 - Use for retesting if results are incorrect
 - Use for QC, validation, or training



Supply requests

- Check the inventory of kits and outdates to reduce the number of shipments.
 - BVBlue, OSOM rapid Trich, pH paper, urine dipsticks, controls
- Send your request 1 month in advance to Wayne Hall, Lorna Rabe, and Edward (Ted) Livant
- Include at least one other person from your site on the e-mail request incase the initial person is not available to receive the shipment.
- MTN does not have supplies in stock. We order from the US distributor. Kits are delivered to our lab in 1-2 weeks depending on the item and if they are in stock at the distributor



Invoices for supplies

- the invoice sent with the package is an estimate of the cost for customs purposes
- A second invoice is sent from the MTN finance department about 1 month later that includes the shipping costs
- Any billing inquiries are to be directed to:
Nicole Lazor: lazorn@mwri.magee.edu

How to store Gram stain slides and enter into LDMS

- Use two slide boxes, one for the slides to be shipped and the other for the backup slide
- Entering into LDMS:
 - storage module
 - Configuration
 - Fill order: change to top to bottom, left to right
 - Coordinates: columns, rows
 - # of columns: 1, # of rows: 100



Shipping slides to Network Lab

- Ship only boxes of 100 slides
- Notify Ted, Wayne, and Lorna when shipping slides
 - Send tracking number
 - LDMS formatted manifest: send by e-mail (not .csv)
- Preparing the slides for shipment
 - Put several paper towels in the lid of the boxes to prevent slides from moving. (shake the box and if it rattles add more paper towels.
 - Wrap tape around the box to secure lid.
 - Use bubble wrap, paper, or Styrofoam packing material
- FedEx, DHL, or World Courier



BVBlue

- Rapid test that detects elevated levels of sialidase
- Previous published studies showed a sensitivity between 89-100%
- Preliminary results from the VOICE enrollment visits

Results: BVBlue vs Gram Stain in VOICE Study

BVBlue test N=371	Gram stain	
	Positive (Nugent 7-10) n = 151	Negative (Nugent 0-6) n = 220
Positive	78	12
Negative	73	208

Sensitivity 52%

Specificity 96%

Positive Predictive value 87%

Negative Predictive value 74%

BVBlue Results by Sites

Sites	Total tested	BVBlue positive	Gram stain score 7-10	Sensitivity	Specificity
All Sites	371	90 (24%)	151 (41%)	52%	96%
Verulam	94	23 (24%)	41 (43%)	49%	94%
Botha Hills	67	18 (27%)	35 (52%)	49%	97%
Umkomaas	59	7 (12%)	14 (24%)	50%	100%
RK Khan	53	11 (21%)	22 (42%)	50%	97%
Tongaat	33	14 (42%)	13 (39%)	69%	75%
Overport	18	3 (17%)	9 (50%)	33%	100%
Isipingo/eThikwini	8	2 (25%)	3 (38%)	67%	100%
Aurum/WRHI	27	7 (26%)	10 (37%)	60%	94%
Zimbabwe	8	1 (13%)	2 (25%)	50%	100%
Uganda	4	3 (75%)	2 (50%)	100%	50%

Women with Gram Stain Positive for BV (Score 7-10)

	BVBlue pos	BVBlue neg	P value
N	78	73	
Mean Gram stain score	8.7	8.3	0.03 ¹
Yeast present on GS	4 (5%)	4 (5%)	NS

	WBC²				P value=0.009 ³
	0	1+	2+	3+	
BVBlue pos	51	17	10	0	
BVBlue neg	37	18	9	9	

1. Mann Whitney

2. <1 WBC/1000X=1+, 1-4 WBC/1000X=2+, 5-30 WBC/1000X=3+

3. Chi Square for linear trend




Diagnosis of BV

- ❑ BVBlue will continue in VOICE
- ❑ BVBlue will NOT be used in MTN 020 (ASPIRE) study
- ❑ Amsel criteria will be use for MTN 020



Wet mount proficiency test


- Send list of readers and their email addresses to Michele Austin and Lorna Rabe
- Please stress to the readers that they need to take the tests as soon as possible.
 - If they fail the initial test they will need to take another test before the deadline date
- A set of training slides are available at the MTN website located in the network/laboratory/wet mount proficiency
 - Use this as a refresher or to supplement your training



Distinguishable
Epithelial Cells (5)
(Nuclei are
present)

Clue Cells (3)

Note: This
field is positive
for clue cells.
 $3/5 = 60\%$





Lactobacillus


Note: This slide is negative for clue cells.
Notice the large number of neutrophils.

Reading KOH Wet Preps

- The following tips can be utilized in order to correctly identify pseudohyphae and yeast in a wet prep.
 1. Pseudohyphae must have yeast “budding off” of the tubular structure in order for it to be correctly identified as pseudohyphae. If there are no “buds” present, it is not pseudohyphae.
 2. In order for yeast not to be mistaken for nuclei of cells, they must be “budding”. This form resembles a larger oval with a smaller oval attached at one end. Without the bud, the yeast would look like a nuclei.


Example:





Artifacts

No yeast budding off of the structures are observed. This can not be called pseudohyphae.



Note: Please note the size of the epithelial cells in relation to the actual pseudohyphae and the artifact.

Pseudohyphae

Epithelial cell

Epithelial cell

Budding yeast

Artifact

Note: No budding yeast are observed coming off of the structure. The size is too large and the shape is not consistent with pseudohyphae.