



UZ-UCSF Laboratory Rapid HIV Testing QC plan for HPTN 035

Yaw Agyei

HPTN LAB QA/QC Coordinator

Purpose

- Ensure that sites are adhering to written procedures when performing rapid HIV testing and that the results are accurate.

Procedure

- All specimens tested with Rapid HIV tests at the site lab, will be sent to UZ-UCSF Central Lab
- UZ-UCSF Central Lab will retest every 10th serum/plasma specimen, depending on the type of sample collected for the particular visit.
- Test used will be FDA approved ELISA (BIO-RAD HIV 1/2 O plus)

Reporting QC checks

- Results from QC checks will be available for review at both the sites and UZ-UCSF central laboratory.
- Where there are discordant results between site and central lab – the rapid HIV algorithm (Oraquick and Determine) will be repeated.
- If still discordant, an ELISA, then a Western Blot test will be performed.

e.g.


Site Results

HIV +ve (both rapid tests +ve)

Central Lab

ELISA HIV -ve

Repeat Oraquick and Determine



If both +ve /discordant



Western Blot

QA Checks

- Site supervisor will perform monthly spot checks (using a check–list) to observe adherence to SOP's and QA procedures.

Check list questions (testing done by Med Lab scientist)

- Testing algorithm used
- Presence of SOP
- Procedure followed
- Expiration date of kit recorded on worksheet
- Test devices labelled correctly
- Kits logged and stored correctly
- Safety procedures and waste disposal procedures being followed

Any Questions or Comments?