



MTN 017 Laboratory Training

Pam Kunjara, MT ASCP
MTN Laboratory Center
Magee-Womens Research Institute
Pittsburgh, PA

Objectives

- ◆ Overview of Screening Lab testing
- ◆ HIV Testing Algorithm
- ◆ HIV Result CRF
- ◆ Q&A






SSP Lab Section

- ◆ This will be your best resource for MTN 017 lab questions
- ◆ Make sure you have the most current version
 - ◆ Available at www.mtnstopshiv.org

MTN Laboratory Contacts

- ◆ Laboratory Center (LC)

- ◆ Pam Kunjara 412-641-6393
pkunjara@mwri.magee.edu



- ◆ JHU Clinical Pharmacology Analytical Lab (JHU CPAL)

- ◆ Craig Hendrix chendrix@jhmi.edu
- ◆ Mark Marzinke mmarzin1@jhmi.edu

Overview of Lab Testing by Visit

	VST 1 SCR	VST 2 ENR	VST 3 MID	VST 4 END	VST 5 PD2	VST 6 MID	VST 7 END	VST 8 PD3	VST 9 MID	VST 10 END
Rectal HSV detection	★	★	★	★	★	★	★	★	★	★
Anal HPV		X								
Rectal GC/CT	X	X	★	X	★	★	X	★	★	X
Rectal sponge for PK		X	X	X	X	X	X	X	X	X
Rectal sponge for PD			X	X	X	X	X	X	X	X
Rectal sponge immuno*		X		X			X			X
Rectal biopsy Proteomics*		X		X			X			X
Rectal biopsy Histology*		X		X			X			X
Rectal biopsies Pheno*		X		X			X			X
Rectal biopsy Gene Array*		X		X			X			X
Rectal biopsies PD*		X		X			X			X
Rectal biopsies PK*				X			X			X

*Tissue subset only



Overview of Lab Testing


Urine Specimens

- ◆ Urinalysis
- ◆ NAAT for GC/CT
 - ◆ GenProbe Aptima
 - ◆ GeneXpert
- ◆ *Separate urine for Urinalysis.
Refrigerate urine for NAAT.*

Overview of Lab Testing

Blood Specimens

- ◆ Red top or Green top tube (Na Hep)
 - ◆ Chemistries
 - ◆ ALT
 - ◆ AST
 - ◆ Creatinine – Calculation tool for Creatinine clearance available www.mtnstopshiv.org
 - ◆ Hepatitis B Surface Antigen and Antibody
 - ◆ Hepatitis C Antibody
 - ◆ Syphilis Serology
 - ◆ HSV Testing





Overview of Lab Testing

Blood Specimens

- 💧 Whole blood: blue top tube (Na Citrate)
 - 💧 PT/INR
- 💧 Whole blood: purple top tube (EDTA)
 - 💧 Complete Blood Count (CBC) with diff and platelets
 - 💧 HIV Testing – follow algorithm (Protocol Appendix II)



HIV Testing Algorithm



*CLIA certified labs may perform 1 rapid test
Ind: Indeterminate test results
EIA: Enzyme Immunoassay

SAMPLE: DO NOT FAX TO DATAFAX
MTN-017 (198)



Visit Code **1**

Participant ID
Site Number Participant Number Ck

Initial Specimen Collection Date
dd MMM yy

HIV Results

1. Rapid HIV test 1:	Not done/Not collected <input type="checkbox"/>	kit code <input type="text"/> <input type="text"/>	Alternate Collection Date dd <input type="text"/> <input type="text"/> MMM <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> yy <input type="text"/> <input type="text"/>	negative <input type="checkbox"/>	positive <input type="checkbox"/>	
2. Rapid HIV test 2:	Not done/Not collected <input type="checkbox"/>	kit code <input type="text"/> <input type="text"/>	Alternate Collection Date dd <input type="text"/> <input type="text"/> MMM <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> yy <input type="text"/> <input type="text"/>	negative <input type="checkbox"/>	positive <input type="checkbox"/>	
3. HIV-EIA:	Not done/Not collected <input type="checkbox"/>		Alternate Collection Date dd <input type="text"/> <input type="text"/> MMM <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> yy <input type="text"/> <input type="text"/>	negative <input type="checkbox"/>	positive <input type="checkbox"/>	indeterminate <input type="checkbox"/>
<p><i>If any are positive at Screening or Enrollment, participant is ineligible. If any are positive during follow-up, complete HIV Confirmatory Results form and Clinical Product Hold/Discontinuation Log.</i></p> <p><i>If indeterminate at Screening, participant is ineligible. If indeterminate during follow-up, complete HIV Confirmatory Results form and Clinical Product Hold/Discontinuation Log.</i></p>						
4. Was plasma stored for the Network Lab (Confirmatory HIV Serology)?	yes <input type="checkbox"/>	no <input type="checkbox"/>	not required <input type="checkbox"/>	If no or not required, end of form.		
<p>Note: Plasma storage for the Network Lab is required at all End Period Visits, any time there is a positive rapid or EIA HIV test result during follow-up, and any time Sample 2 is collected.</p>						
4a. Plasma for storage collection date:			dd <input type="text"/> <input type="text"/>	MMM <input type="text"/> <input type="text"/> <input type="text"/>	yy <input type="text"/> <input type="text"/>	

Comments:

HIV Results (HIV-1)

Purpose: This form is used to document the participant's HIV rapid test or EIA results, and plasma storage for HIV confirmatory testing as specified in the protocol.

General Information/Instructions: Record test results on this form as they become available. Fax this form into SCHARP DataFax once results for all collected specimens are recorded on the form.

- Initial Specimen Collection Date:** Record the date that the first specimen was collected (NOT the date results were reported or recorded on the form). A complete date is required.
- Alternate Collection Date:** This date is to be completed ONLY if the specimen was collected on a date after the Initial Specimen Collection Date. A specimen collected for the same visit but on a different date should be recorded on the same form.
- Not done/Not collected:** Mark this box in the event that a specimen was not collected or if the specimen was collected, but a result is not available due to specimen loss or damage. Record the reason why the result is not available in the Comments.

Item-specific Instructions:

Visit Code: Record the visit code assigned to this visit. Refer to the Study-Specific Procedures (SSP) Manual for more specific information on assigning visit codes.

Items 1 and 2: Record the assigned two-digit rapid test kit code. *Note: More test kit codes may be added to the list as the study proceeds.*

Kit Code	Rapid Test
01	Determine
02	OraQuick
03	Uni-Gold Recombigen
04	Bioline
05	Clearview Statpak

If item 1, 2, or 3 is positive (meaning the participant had at least one positive HIV test result) during follow-up, complete a new HIV Confirmatory Results form and a Clinical Product Hold/Discontinuation (PH) Log.

Item 3: If positive during follow-up, complete a new HIV Confirmatory Results (HCR) form and a Clinical Product Hold/Discontinuation (PH) Log.



Collection of Rectal Specimens

- ◆ NAAT for GC/CT
 - ◆ *Product gel may cause interference during testing. Please be sure to avoid gel during collection.*
 - ◆ Use Gene Xpert collection swab.
 - ◆ After collection place the swab in the transport tube, break off shaft and cap. The specimen can now remain at 2-30°C for 30 days.
 - ◆ The results are sent to the clinic and are reported on a CRF.

Any Questions?

