

Assay	Turn Around Time Requirements
Urine pregnancy test	A
Urine SDA for gonorrhea and Chlamydia (neat method)	B
Dipstick Urinalysis	A
HIV Rapid Tests	A
HIV Western Blot	B
Complete blood count	C
Serum Chemistries (AST, ALT, Creatinine, Phosphate)	C
Hepatitis B Surface Antigen and Antibody	B
Syphilis Serology	B
CD4+ T Cell Count	C
HIV-1 RNA PCR	B
Plasma Archive	D
Pap Smear	B
Vaginal pH	E
KOH wet preparation	C
Trichomonas Rapid Test	A
BV Rapid Test	A
Vaginal Gram Stain	Slides made at bedside, batched and ship per NL instructions
Vaginal/Endocervical Swabs	Freeze within 8 Hours

Category A: Must be done during visit when collected because may affect eligibility requirements or cause product hold

Category B: Batched at site discretion and within site SOP's for specimen stability-results required for enrollment and/or clinical management

Category C: Per site SOP for specimen stability-generally same day or less than a week.

Category D: On ice or refrigerated: transport to lab and freeze within 8 hours of collection. Room temperature: transport to lab and freeze within 4 hours of collection.

Store at site until notified by Study Team or Network Lab.

Category E: Point of Care