

AMPLIPREP-AMPLICOR VIRAL PCR READINGS VS TAQMAN VIRAL PCR READINGS

NO.	Amplicor reading 1 Copies/ml	Log value	Amplicor reading 2 Copies/ml	Log value	TaqMan reading 1 Copies/ml	Log value	TaqMan reading 2 copies/ml	Log value
1	1060 000	6.03	1420 000	6.15	698 162	5.84	523 020	5.71
2	45700	4.66	59800	4.78	236980	5.37	328810	5.51
3	451000	5.65	453000	5.66	458520	5.66	458520	5.66
4	76600	4.88	175000	5.24	110476	5.04	110476	5.04
5	521000	5.72	554000	5.74	818560	5.91	818560	5.91
6	288000	5.46	241000	5.38	35537	4.55	192370	5.28
7	60500	4.78	49000	4.69	18087	4.26	11750	4.07
8	390000	5.59	709000	5.85	640264	5.81	773524	5.89
9	20600	4.31	28700	4.46	8903	3.95	11410	4.06
10	8930	3.95	25100	4.39	51403	4.71	46595	4.67
11	4280	3.63	7590	3.88	20782	4.32	29450	4.47
12	1460	3.16	657	2.82	980	2.99	1130	3.05
13	34300	4.53	33700	4.53	28847	4.46	34270	4.54
14	214000	5.33	573000	5.76	250543	5.39	155660	5.19
15	186000	5.27	68300	4.83	475573	5.67	475573	5.67
16	193000	5.28	433000	5.64	271538	5.43	280950	5.45
17	256000	5.4	208000	5.32	71476	4.85	71476	4.85
18	<400		<400		<40		<40	
19	15600	4.19	10600	4.02	13068	4.11	13670	4.14
20	211000	5.32	89900	4.95	98567	4.99	98567	4.99
21	255000	5.41	277000	5.44	203088	5.31	176620	5.25
22	267000	5.43	271000	5.43	260000	5.41	260000	5.41
23	726000	5.86	650000	5.81	376533	5.58	376533	5.58
24	10800	4.03	94900	4.98	90000	4.95	90000	4.95
25	171000	5.23	350000	5.54	267942	5.43	267942	5.43
26	360000	5.56	301000	5.48	352875	5.55	352875	5.55
27	33000	4.52	27500	4.44	18456	4.27	15580	4.19
28	300000	5.48	138000	5.14	82555	4.91	92585	4.97
29	599000	5.78	524000	5.72	252857	5.4	152110	5.18
30	773000	5.89	750000	5.88	494953	5.69	498119	5.69
31	<400		<400		<40		<40	
32	<400		<400		<40		<40	
33	6700	3.83	6700	3.83	12495	4.09	12495	4.09
34	223000	5.35	223000	5.35	76695	4.88	76695	4.88
35	92000	4.96	85400	4.93	81125	4.91	81125	4.91
36	64800	4.81	33900	4.53	15843	4.19	23588	4.37
37	<400		<400		<40		<40	
38	<400		<400		94	1.97	97	1.99
39	<400		<400		<40		<40	
40	<400		<400		<40		<40	

PROTOCOL SAMPLES RUN ON AMPLIPREP-AMPLICOR VS TAQMAN

		AMPLIPREP- AMPLICOR Copies/ml	LOG VALUE	TAQMAN Copies/ml	LOG VALUE
1	53120192	1060000	6.03	698162	5.84

2	53121202	45700	4.66	236980	5.37
3	53121320	451000	5.65	453000	5.66
4	53121368	175000	5.24	110476	5.04
5	3121159	521000	5.72	818560	5.91
6	3121160	288000	5.46	192370	5.28
7	3121161	60500	4.78	18987	4.28
8	53121386	390000	5.59	640264	5.8
9	3121168	20600	4.31	11410	4.05
10	3121163	25100	4.39	46595	4.67